





The Creation Of The Calibration Methods And The Analysis Of Alcoholic Drinks Using 3-Pentanol Or Ethanol As Internal Standards (ISTD) in Agilent ChemStation B.04.03



Minsk, 2025

1. The creation of calibration methods

- 1.1. Acquire the chromatograms of standard solution for calibration (SScalibr).
- 1.2. Click “Method” in the main menu and press “Load Method...” in the list.
- 1.3. Select any method and then save the method with new name: “Traditional ISTD” in the case of the using 3-pentanol as internal standard and “Ethanol ISTD” in the case of the using ethanol as internal standard.
- 1.4. Set and integrate the peaks of the volatiles compounds in the chromatograms.



$$C[\text{mg/L AA}] = C[\mu\text{g/g}] \cdot \rho[\text{g/L}] / (10 \cdot ABV)$$

$$C_{\text{ethanol}} = \rho_{\text{ethanol}} = 789270 \text{ mg/L AA}$$

1. The creation of calibration methods

- 1.5.** Fill the input lines in dialog window “Calibration Setting”. Fill in the cell “Amount Units”: “μg/g” in the case of “Traditional ISTD” project and “mg/L AA” in the case of “Ethanol ISTD” project.
- 1.6.** Input the names of the compounds in “Compound” column and the certified values of compound concentration in “Amt[μg/g]” or “Amt[mg/L AA]” column. It will be named these concentration values of compound i (i = 1, ..., 13) in SScalibr as (μg/g) in the case of 3-pentanol as ISTD and (mg/L AA) in the case of ethanol as ISTD. Concentration of ethanol in mg/L AA is equal to ethanol density (= 789270 mg/L).

Calibration Report Batch View Abt

New Calibration Table...
Delete Calibration Table...
Recalibrate...
Add Level...
Add Peaks...
Calibration Settings...

Title: Calibration 3-Pentanol ISTD

Default RT Windows
Minutes %
Reference Peaks 0.00 + 5.00
Other Peaks 0.00 5.00

Default Calibration Curve
Type Linear
Origin Force
Weight Equal

Amount Units: μg/g


Title: Calibration Ethanol ISTD

Default RT Windows
Minutes %
Reference Peaks 0.00 + 5.00
Other Peaks 0.00 5.00

Default Calibration Curve
Type Linear
Origin Force
Weight Equal

Amount Units: mg/L AA

Parameters



Calibration Table

#	RT	Compound	ISTD	Lvl	Amt[μg/g]
1	6.2	acetaldehyde	No	1	180.68
2	8.6	methyl acetate	No	1	287.96
3	10.5	ethyl acetate	No	1	260.11
4	10.7	acetal	No	1	233.31
5	11.6	methanol	No	1	203.11
6	13.9	ethanol	No	1	333000.00
7	18.4	2-butanol	No	1	210.07
8	19.3	n-propanol	No	1	204.23
9	23.1	2-methyl-1-propanol	No	1	200.14
10	24.2	3-pentanol	Yes	1	226.16
11	26.9	n-butanol	No	1	209.60
12	31.2	2-methyl-1-butanol	Yes	1	216.87
13	31.3	3-methyl-1-butanol	No	1	207.88

Calibration tables

Calibration Table

#	RT	Signal	Compound	ISTD	Lvl	Amt[mg/L AA]
1	6.181	FID1 A	acetaldehyde	No	1	428.250
2	8.554	FID1 A	methyl acetate	No	1	682.510
3	10.471	FID1 A	ethyl acetate	No	1	616.510
4	10.693	FID1 A	acetal	No	1	552.980
5	11.557	FID1 A	methanol	No	1	481.410
6	13.912	FID1 A	ethanol	Yes	1	789270.000
7	18.434	FID1 A	2-butanol	No	1	497.890
8	19.318	FID1 A	n-propanol	Yes	1	484.060
9	23.118	FID1 A	2-methyl-1-propanol	No	1	474.370
10	24.168	FID1 A	3-pentanol	No	1	536.030
11	26.872	FID1 A	n-butanol	No	1	496.770
12	31.155	FID1 A	2-methyl-1-butanol	No	1	514.010
13	31.348	FID1 A	3-methyl-1-butanol	No	1	492.710

$$C[\text{mg/L AA}] = C[\mu\text{g/g}] \cdot \rho[\text{g/L}] / (10 \cdot \text{ABV})$$

$$C_{\text{ethanol}} = \rho_{\text{ethanol}} = 789270 \text{ mg/L AA}$$

The calibration tables in the methods: “Traditional ISTD” (3-pentanol as ISTD) and “Ethanol ISTD”

#	RT	Compound	ISTD	Implem.#	Amt[µg/g]	Area	Rsp.Factor
1	6.2	acetaldehyde	No	1	180.68	2.82e-4	640847
				2	180.68	2.86e-4	630730
				3	180.68	2.84e-4	637261
2	8.6	methyl acetate	No	1	287.96	3.67e-4	785668
				2	287.96	3.69e-4	779607
				3	287.96	3.62e-4	795822
3	10.5	ethyl acetate	No	1	260.11	4.64e-4	561140
				2	260.11	4.64e-4	560645
				3	260.11	4.57e-4	568818
4	10.7	acetal	No	1	233.31	5.56e-4	419903
				2	233.31	5.60e-4	416793
				3	233.31	5.49e-4	425099
5	11.6	methanol	No	1	203.11	3.37e-4	602347
				2	203.11	3.39e-4	599041
				3	203.11	3.35e-4	606212
6	13.9	ethanol	No	2	333000.00	0.649	513446
				3	333000.00	0.639	520885
				1	333000.00	0.644	517252
7	18.4	2-butanol	No	1	210.07	6.40e-4	328216
				2	210.07	6.48e-4	324187
				3	210.07	6.38e-4	329344
8	19.3	n-propanol	No	1	204.23	6.11e-4	334346
				2	204.23	6.14e-4	332515
				3	204.23	6.03e-4	338490
9	23.1	2-methyl-1-propanol	No	1	200.14	7.06e-4	283679
				2	200.14	7.16e-4	279593
				3	200.14	7.02e-4	285178
10	24.2	3-pentanol	Yes	1	226.16	7.60e-4	297519
				2	226.16	7.68e-4	294411
				3	226.16	7.54e-4	299922
11	26.9	n-butanol	No	1	209.60	6.87e-4	304936
				2	209.60	6.93e-4	302559
				3	209.60	6.83e-4	306882
12	31.2	2-methyl-1-butanol	No	1	216.87	7.71e-4	281380
				2	216.87	7.75e-4	279748
				3	216.87	7.67e-4	282590
13	31.3	3-methyl-1-butanol	No	1	207.88	7.30e-4	284585
				2	207.88	7.32e-4	283981
				3	207.88	7.25e-4	286813

#	RT	Compound	ISTD	Implem.#	Amt[mg/L AA]	Area	Rsp.Factor
1	6.2	acetaldehyde	No	1	428.25	2.82e-4	1518944
				2	428.25	2.86e-4	1494965
				3	428.25	2.84e-4	1510443
2	8.5	methyl acetate	No	1	682.51	3.67e-4	1862156
				2	682.51	3.69e-4	1847790
				3	682.51	3.62e-4	1886223
3	10.5	ethyl acetate	No	1	616.51	4.64e-4	1329996
				2	616.51	4.64e-4	1328834
				3	616.51	4.57e-4	1348206
4	10.7	acetal	No	1	552.98	5.56e-4	995234
				2	552.98	5.60e-4	987862
				3	552.98	5.49e-4	1007548
5	11.6	methanol	No	1	481.41	3.37e-4	1427679
				2	481.41	3.39e-4	1419844
				3	481.41	3.35e-4	1436840
6	13.9	ethanol	Yes	1	789270.00	0.644	1225966
				2	789270.00	0.649	1216959
				3	789270.00	0.639	1234591
7	18.4	2-butanol	No	1	497.89	6.40e-4	777910
				2	497.89	6.48e-4	768361
				3	497.89	6.38e-4	780582
8	19.3	n-propanol	No	1	484.06	6.11e-4	792457
				2	484.06	6.14e-4	788117
				3	484.06	6.03e-4	802278
9	23.1	2-methyl-1-propanol	No	1	474.37	7.06e-4	672372
				2	474.37	7.16e-4	662689
				3	474.37	7.02e-4	675926
10	24.2	3-pentanol	No	1	536.03	7.60e-4	705161
				2	536.03	7.68e-4	697794
				3	536.03	7.54e-4	710856
11	26.9	n-butanol	No	1	496.77	6.87e-4	722724
				2	496.77	6.93e-4	717091
				3	496.77	6.83e-4	727336
12	31.2	2-methyl-1-butanol	No	1	514.01	7.71e-4	666907
				2	514.01	7.75e-4	663039
				3	514.01	7.67e-4	669775
13	31.3	3-methyl-1-butanol	No	1	492.71	7.30e-4	674513
				2	492.71	7.32e-4	673081
				3	492.71	7.25e-4	679794


2. Analysis of alcoholic beverages

- 2.1. Acquire the chromatograms of samples of alcoholic beverages. Peaks of compounds are identified and integrated, including the ethanol peaks.
- 2.2. Click “Report” in the main menu and press “Specify Report...”. Set in “Quantitative Results”: “Calculate” – “ISTDed on” – “Area”.


Samples of alcoholic beverages



Chromatograms



Specify Report



Report (traditional internal standard)

Sorted By	:	Signal
Calib. Data Modified	:	Thursday, December 26, 2024 7:15:13 PM
Multiplier	:	1.0000
Dilution	:	1.0000
Use Multiplier & Dilution Factor with ISTDs		
Sample ISTD Information:		
ISTD #	ISTD Amount [µg/g]	Name
1	227.80000	3-pentanol

Signal 1: FID1 A, Front Signal

RetTime [min]	Type	ISTD used	Area [V*s]	Amt/Area ratio	Amount [µg/g]	Grp	Name
6.185	BB	1	9.02712e-5	2.14029	58.17736		acetaldehyde
8.534	BB	1	1.48048e-5	2.64738	11.80194		methyl acetate
10.446	BB	1	8.49700e-4	1.89561	485.00482		ethyl acetate
10.668	BB	1	1.39654e-4	1.41479	59.49469		acetal
11.523	BB	1	4.37130e-3	2.02682	2667.82824		methanol
13.930	BB	1	7.71661e-1	1.73975	4.04246e5		ethanol
18.434		1	-	-	-		2-butanol
19.271	BB	1	4.53447e-4	1.12726	153.91588		n-propanol
23.052	BB	1	3.04911e-4	9.51327e-1	87.34440		2-methyl-1-propanol
24.102	BB	I	7.56521e-4	1.00000	227.80000		3-pentanol
26.802	BB	1	9.14371e-5	1.02527	28.22876		n-butanol
31.091	BB	1	1.38481e-4	9.46042e-1	39.44860		2-methyl-1-butanol
31.284	BB	1	6.80836e-4	9.59114e-1	196.62790		3-methyl-1-butanol

$$C^{\text{measured}} [\mu\text{g/g}] = C_{3\text{-pentanol}} \cdot RRF^{3\text{-pentanol}} \cdot \frac{A}{A_{3\text{-pentanol}}}$$

Report (ethanol as internal standard)

Sorted By	:	Signal
Calib. Data Modified	:	Tuesday, December 24, 2024 7:10:13 PM
Multiplier	:	1.0000
Dilution	:	1.0000
Use Multiplier & Dilution Factor with ISTDs		
Sample ISTD Information:		
ISTD #	ISTD Amount [mg/L AA]	Name
1	7.89270e5	ethanol

Signal 1: FID1 A, Front Signal

RetTime [min]	Type	ISTD used	Area [V*s]	Amt/Area ratio	Amount [mg/L AA]	Grp	Name
6.185	BB	1	9.02712e-5	1.23025	113.59049		acetaldehyde
8.534	BB	1	1.48048e-5	1.52169	23.04247		methyl acetate
10.446	BB	1	8.49700e-4	1.08959	946.95446		ethyl acetate
10.668	BB	1	1.39654e-4	8.13209e-1	116.15915		acetal
11.523	BB	1	4.37130e-3	1.16502	5208.86353		methanol
13.930	BB	I	7.71661e-1	1.00000	7.89270e5		ethanol
18.431		1	-	-	-		2-butanol
19.271	BB	1	4.53447e-4	6.47943e-1	300.51277		n-propanol
23.052	BB	1	3.04911e-4	5.46821e-1	170.53630		2-methyl-1-propanol
24.102	BB	1	7.56521e-4	5.74786e-1	444.76045		3-pentanol
26.802	BB	1	9.14371e-5	5.89298e-1	55.11329		n-butanol
31.091	BB	1	1.38481e-4	5.43773e-1	77.02034		2-methyl-1-butanol
31.284	BB	1	6.80836e-4	5.51296e-1	383.90693		3-methyl-1-butanol

$$C^{\text{measured}} [\text{mg/L AA}] = \rho_{\text{ethanol}} \cdot RRF^{\text{ethanol}} \cdot \frac{A}{A_{\text{ethanol}}}$$

Single-level calibration with internal standard – ethanol

1. Certified values of concentration in mg/L AA (per anhydrous alcohol)

$$C[\text{mg/L AA}] = \frac{C[\mu\text{g/g}] \cdot \rho[\text{g/L}]}{10 \cdot \text{ABV}}$$


One calibration solution (SS-C) for RRF^{ethanol} determination



Samples of alcoholic beverages



2. Calculation of the relative response factor RRF^{ethanol}



Calibration with SS-C using ethanol as internal standard

$$RRF^{\text{ethanol}} = \frac{C^{\text{certified}}}{\rho_{\text{ethanol}}} / \frac{A}{A_{\text{ethanol}}}$$

$$\rho_{\text{ethanol}} = 789270 \text{ mg/L}$$

no need to add an internal standard substance

3. Calculation of compound concentration in the products

$$C^{\text{measured}} [\text{mg/L AA}] = \rho_{\text{ethanol}} \cdot RRF^{\text{ethanol}} \cdot \frac{A}{A_{\text{ethanol}}}$$

The concentration of volatile compounds is calculated in mg/L AA
There is **no need to determine the strength** of the alcoholic beverage samples.