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\*Charapitsa S., et al. Direct Determination of Volatile Compounds in Spirit Drinks by Gas Chromatography.

*Journal of Agricultural and Food Chemistry*, 2013, 61, 2950-2956. [doi.org/10.1021/jf3044956](https://doi.org/10.1021/jf3044956)

Charapitsa S., et al. The establishment of metrological characteristics of the method “Ethanol as Internal Standard” for the direct determination of volatile compounds in alcoholic products .

*Journal of Chemical Metrology*, 2018, 12, 59-69. [doi.org/10.25135/jcm.14.18.02.063](https://doi.org/10.25135/jcm.14.18.02.063)

Charapitsa S., et al Single-Laboratory Validation of a Gas Chromatographic Method of Direct Determination of Volatile Compounds in Spirit Drinks:

Need for an Improved Interlaboratory Study, *Journal of AOAC International*, 2019, 102, 2, 669-672. [doi.org/10.5740/jaoacint.18-0258](https://doi.org/10.5740/jaoacint.18-0258)

Charapitsa S., et al. Interlaboratory study of ethanol usage as an internal standard in direct determination of volatile compounds in alcoholic products,

*42nd World Congress of Vine and Wine*, 2019, [doi.org/10.1051/bioconf/20191502030](https://doi.org/10.1051/bioconf/20191502030)

Charapitsa S., et al. . The Perspectives of Ethanol Usage as an Internal Standard for the Quantification of Volatile Compounds in Alcoholic Products by GC-MS.

*Journal of Mass Spectrometry*, 2019, [doi.org/10.1002/jms.4493](https://doi.org/10.1002/jms.4493)

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The presence of these **three additional terms (4, 5, 6)** indicates a higher uncertainty during the determination of concentrations of volatile compounds by the traditional method than “Ethanol is Internal Standard” method.

These additional terms are responsible for the uncertainty in the determination of concentration of internal standard, density and ABV value of a test sample.

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## SHORT COMMUNICATION

# Single-Laboratory Validation of a Gas Chromatographic Method of Direct Determination of Volatile Compounds in Spirit Drinks: **Need for an Improved Interlaboratory Study**

SIARHEI V. CHARAPITSA and SVETLANA N. SYTOVA

Belarusian State University, Institute for Nuclear Problems, Bobruyskaya St, 11, Minsk, Belarus

ANTON L. KORBAN and LIDIA N. SOBOLENKO

Belarusian State University, Institute for Nuclear Problems, Bobruyskaya St, 11, Minsk, Belarus; Belarusian State University, Nezavisimosti Ave, 4, Minsk, Belarus

**Background:** The quality and safety control of an alcoholic drink is mainly the establishment of its chemical content, particularly the quantity of volatile compounds. **Objective:** A single-laboratory validation of a gas chromatographic method of direct determination of volatile compounds in spirit drinks was conducted. The discussed method applies ethanol, the major volatile component of an alcoholic beverage, as an internal standard. Possible algorithms of method validation based on interlaboratory study were proposed and described. **Methods:** Seven standard solutions of the following volatile compounds were prepared gravimetrically in 40% (v/v) water–ethanol solution: acetaldehyde, methyl acetate, ethyl acetate, methanol, 2-propanol, 1-propanol, isobutanol, 1-butanol, and isoamylol. Each sample was measured with the proposed method 30 times in repeatability conditions. **Results:** Flame ionization detector response was linearly

in micrograms per gram units (1). It is important to note that legislative documents establish final presentation of volatile compound concentrations in milligrams per liter absolute alcohol (AA; or grams/hectoliter AA) units (2–4). To present concentrations in these quantity dimensions, one should perform an additional procedure of measuring the alcohol by volume content in the test sample.

## Methods


The current study is a further exploration of our previously published “Ethanol as Internal Standard” method (5–6). It consists of the determination of Relative Response Factor (RRF) for ethanol. First, response factor for *i*-th volatile and ethanol are determined according to the following formula:

$$\text{Response Factor}_i = \frac{\text{Concn}_{i/eth}^{st}}{\text{Area}_{i/eth}^{st}}, \quad (1)$$

where  $\text{Concn}_i^{st}$  is the concentration of *i*-th volatile or ethanol

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Data File C:\CHEM32\1\DATA\ALCOHOL\BRANDY.D  
Sample Name: RVE-K-1

```

=====
Acq. Operator   : shest                      Seq. Line : 21
Acq. Instrument : Instrument 1                Location  : Vial 106
Injection Date  : 15.02.2019 0:26:17         Inj       : 1
                                           Inj Volume: 1 µl

Acq. Method    : C:\HPCHEM\1\METHODS\1\SPIRT_RV.M
Last changed   : 14.02.2019 21:05:08 by shest
Analysis Method : C:\CHEM32\1\METHODS\EXTERNAL_STANDARD_RV.M
Last changed   : 13.05.2019 13:53:13
                (modified after loading)
Additional Info : Peak(s) manually integrated
=====
External Standard Report
=====
Sorted By      : Signal
Calib. Data Modified : 13.05.2019 13:52:10
Multiplier    : 1.0000
Dilution      : 1.0000
Use Multiplier & Dilution Factor with ISTDs

Signal 1: FID1 A,

```

RetTime [min]	Type	Area [pA*s]	Amt/Area	Amount [mg/L]	Grp	Name
3.146	BB	19.30616	3.54975	68.53204		Acetaldehyde
3.839	VB	3.89059	3.26799	12.71441		Methyl acetate
4.491	BV	56.77414	1.98393	112.63583		Ethyl acetate
4.651	VB	52.10571	3.33915e-4	1.73989e-2		Methanol, %
5.084		-	-	-		2-propanol
5.297	BB S	1.89385e5	2.09573e-4	39.68984		Ethanol, %
7.770	BB	79.04855	1.36267	107.71724		1-propanol
9.153	BB	422.61612	1.14194	482.60232		2-methyl-1-propanol
10.439	BB	1.11518	1.16797	1.30250		1-butanol
11.792	BB	1213.55713	1.08046	1311.19524		3-methyl-1-butanol
Totals :				2136.40682		

2 Warnings or Errors :

Warning : Calibration warnings (see calibration table listing)  
Warning : Calibrated compound(s) not found

```

=====
*** End of Report ***
=====
Instrument 1 13.05.2019 14:03:45 Page 1 of 1

```

Data File C:\CHEM32\1\DATA\ALCOHOL\BRANDY.D  
Sample Name: RVE-K-1

```

=====
Acq. Operator   : shest                      Seq. Line : 21
Acq. Instrument : Instrument 1                Location  : Vial 106
Injection Date  : 15.02.2019 0:26:17         Inj       : 1
                                           Inj Volume: 1 µl

Acq. Method    : C:\HPCHEM\1\METHODS\1\SPIRT_RV.M
Last changed   : 14.02.2019 21:05:08 by shest
Analysis Method : C:\CHEM32\1\METHODS\ETHANOL_AS_INTERNAL_STANDARD.M
Last changed   : 13.05.2019 14:28:42
                (modified after loading)
Additional Info : Peak(s) manually integrated
=====
Internal Standard Report
=====
Sorted By      : Signal
Calib. Data Modified : 13 May 2019 r. 14:28:42
Multiplier    : 1.0000
Dilution      : 1.0000
Use Multiplier & Dilution Factor with ISTDs
Sample ISTD Information:
ISTD ISTD Amount Name
# [ng/L AA]
-----|-----|-----|-----|-----|-----|-----|
1 7.89300e5 Ethanol

```

Signal 1: FID1 A,

RetTime [min]	Type	ISTD used	Area [pA*s]	Amt/Area ratio	Amount [ng/L AA]	Grp	Name
3.146	BB	1	19.29619	2.13828	171.96268		Acetaldehyde
3.839	VB	1	3.88879	1.96794	31.89511		Methyl acetate
4.491	BV	1	56.76944	1.19541	282.83414		Ethyl acetate
4.651	VB	1	52.09607	1.58959	345.13520		Methanol
5.085		1	-	-	-		2-propanol
5.297	VB S I	1	1.89384e5	1.00000	7.89300e5		Ethanol
7.770	BB	1	79.03934	8.23564e-1	271.29342		1-propanol
9.153	BB	1	422.60141	6.88636e-1	1212.88444		2-methyl-1-propanol
10.439	BB	1	1.11444	7.07093e-1	3.28421		1-butanol
11.792	BB	1	1213.55347	6.52135e-1	3298.33836		3-methyl-1-butanol
Totals without ISTD(s) :					5617.62759		

2 Warnings or Errors :

Warning : Calibration warnings (see calibration table listing)  
Warning : Calibrated compound(s) not found

```

=====
*** End of Report ***
=====
Instrument 1 13.05.2019 14:30:12 Page 1 of 1

```

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42<sup>nd</sup> World Congress of Vine and Wine

<https://doi.org/10.1051/bioconf/20191502030>

## Interlaboratory study of ethanol usage as an internal standard in direct determination of volatile compounds in alcoholic products

S. Charapitsa<sup>1</sup>, S. Sytova<sup>1</sup>, A. Korban<sup>1,2</sup>, L. Sobolenko<sup>1,2</sup>, V. Egorov<sup>2</sup>, S. Leschev<sup>2</sup>, M. Zakharov<sup>3</sup>, R. Čabala<sup>4</sup>, R. Busarova<sup>5</sup>, I. Shestakovich<sup>6</sup>, A. Tolstouhova<sup>6</sup>, S. Ondroušek<sup>7</sup>, J. Vávra<sup>7</sup>, M. Yilmaztekin<sup>8</sup>, and T. Cabaroglu<sup>9</sup>

<sup>1</sup> Institute for Nuclear Problems of Belarusian State University, 220030, Bobruyskaya Str., 11, Minsk, Belarus

<sup>2</sup> Chemistry Faculty, Department of Analytical Chemistry, Belarusian State University, Leningradskaya Str., 14, 220050, Minsk, Belarus

<sup>3</sup> All-Russian Scientific Research Institute of the Brewing, Non-Alcoholic and Wine Industry – Branch of the V.M. Gorbатов Federal Scientific Center of Food Systems of RAS, 119021, Rossolimo Str., 17, Moscow, Russia

<sup>4</sup> Charles University, Faculty of Science, Department of Analytical Chemistry, Hlavova Str., 2030/8, 128 40, Prague 2, Czech Republic

<sup>5</sup> JS “Mosazervinzavod”, 115088, Yuzhnoportovaya Str., 30, Moscow, Russia

<sup>6</sup> Republican Centre for Hygiene, Epidemiology and Public Health, 220099, Kazintsya Str., 50, Minsk, Belarus

<sup>7</sup> General Directorate of Customs, Customs Technical Laboratory, Budějovická Str., 7, 140 96, Prague 4, Czech Republic

<sup>8</sup> Inonu University, Faculty of Engineering, Department of Food Engineering, 44280, Malatya, Turkey

<sup>9</sup> Cukurova University, Faculty of Agriculture, Department of Food Engineering, 01330, Adana, Turkey

**Abstract.** A collaborative interlaboratory study on the method of direct quantitation of volatile compounds in spirit drinks and alcoholic products was conducted. The discussed method applies ethanol, the major volatile component of an alcoholic product, as an internal standard. In this study 9 laboratories from 4 different countries were supplied with standard solutions for gas chromatographic measurements. Five aqueous ethanol 40% (v/v) standard solutions containing target compounds in concentrations ranging from 10 mg/L to 400 mg/L of absolute alcohol were prepared and sent to the participants for quantification of acetaldehyde, methyl acetate, ethyl acetate, methanol, 2-propanol, 1-propanol, 2-methyl-1-propanol, 1-butanol and 3-methyl-1-butanol. The interlaboratory study was evaluated according to the ISO 5725 standards and the Eurachem guide. The within-laboratory precision varied between 0.4% and 7.5% for all samples and compounds, showing a sufficiently high repeatability of the method. The between-laboratory precision was found to vary within a satisfactory range of 0.5% ÷ 10.0%. Precision of the method was well within the range predicted by the Horwitz equation for all analytes. The analysis of trueness showed that the bias of the method is insignificant at the significance level  $\alpha = 5\%$ . The determined concentrations of the analytes compared well to the gravimetric values thus showing very satisfactory accuracy of the method. The results of the interlaboratory study confirmed that “Ethanol as Internal Standard” method is robust and reliable and can be used as a standard reference method for analysing volatile compounds in water-ethanol samples. The possibilities of method validation according to the previously obtained experimental data were shown.

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